

EasyGlide

System**RoMedic**



Instructions for use – English



REF

5050
5020
5060
5070



SWL: 135 - 150 kg/298 - 331 lbs
SWL: 135 - 170 kg/298 - 375 lbs (5050)

EasyGlide

SystemRoMedic™

Instructions for use - English

The sliding board EasyGlide is an ingenious invention. The board helps transfer sitting users to and from wheelchairs, beds, toilets, car seats and many other seated positions. It is especially suitable for levelling out height differences and for bridging short distances.

The rounded corners and the smooth, low-friction top surface provide for easy sliding transfers and make the sliding board very easy to position and remove. The double anti-slip tape on the underside prevents the board from sliding against the underlying surface during transfer.


Check Safety

Visual inspection

Check the condition and function of the assistive device regularly. Verify that the materials are not worn or faded. If there are signs of wear, the product must be discarded.

Stability test

Carefully bend the EasyGlide to check strength and durability. The board should be both elastic and rigid. Discard the EasyGlide when its strength and durability are diminished.



Always read the user manual

Always read the user manuals for all assistive devices used during a transfer.




Keep the user manual where it is accessible to users of the product.

Always make sure that you have the right version of the user manual.


The most recent editions of manuals are available for downloading from our website, www.directhealthcaregroup.com.

Important information

- Do not leave the user unattended during a transfer situation.



Placement/removal of EasyGlide




Place the wheelchair at an angle towards the side of the bed, and as close to the bed as possible. Remove the armrest that is closest to the bed and remove leg supports. Place the EasyGlide like a bridge between the bed and the wheelchair. If it is an EasyGlide model on which one long side can be bent, this should be placed against the wheelchair wheel.

To position the EasyGlide under the client, lean the client slightly in the opposite direction.

When the transfer is complete, EasyGlide is removed by bending upwards.

Using EasyGlide



Active users can transfer themselves with the help of EasyGlide. Grasp the armrest and lean over towards the wheelchair. It is important to remember to lean once more towards the bed when the transfer is complete, so as to be able to move back in the middle of the chair and slide all the way back.

If the client requires support and assistance during the transfer, the EasyGlide should be complemented with an EasyBelt/FlexiBelt.

The caregiver can work from behind or in front of the client.

Tips! Our FootStool can be used to advantage to raise the knees

EasyGlide can also be used for transfers into and out of a car.

Material

Polyethen

Care of the product

The surface of the product can be cleaned with warm water or a soap solution. If disinfection is needed, 70 % ethanol, 45 % isopropanol or similar should be used.

Contact your local distributor if you have any questions about the product and its use.

See www.directhealthcaregroup.com for a complete list of distributors.









Expected lifetime

Up to 10 years of normal use.

Product information

Ref. no.	Size
5020	600 x 330 mm / 23,6 x 12,9"
5050	750 x 330 mm / 29,5 x 12,9"
5060	600 x 250 mm / 23,6 x 9,8"
5070	450 x 200 mm / 17,7 x 7,8"

Symbols

	This product complies with the requirements of the Medical Device Regulation 2017/745		Visual Inspection
	Medical Device		Read the manual
	Caution		Product Code
	Manufacturer information		Batch Code

Simple solutions for great results

SystemRoMedic[®] is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic[®] is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic[®] is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic[®] offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;
www.directhealthcaregroup.com

For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic[®] representative. A complete list of all our partners with their contact details can be found on our website; www.directhealthcaregroup.com.



Moving Health Forward



Medical Device Class I. The product complies with the requirements of the Medical Devices Regulation 2017/745.

Direct Healthcare Group Sverige AB
Torshamnsgatan 35,
SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200
info@directhealthcaregroup.com
www.directhealthcaregroup.com