

INVACARE Action[®] Ampla[™]



The Action Ampla has been specifically developed to take into consideration the body shapes, activity levels and weight limit requirements of Plus-Size users, earning it multiple prestigious design awards. The Ampla has inbuilt transportation bearings, allowing easy movement through doorways or narrow spaces.

The ergonomic backrest allows for more freedom of movement with no obstructions, maximising access to the handrims.

The height adjustable transfer levers can be turned inward or outward to provide support for entering or leaving the chair and can even be lowered to allow table access.

Designed to accommodate the anatomy of both “apple” and “pear” body shapes of Plus Size users, offering comfort and support.

Height adjustable transfer levers.

Multiple active rear wheel positions.

Elastic, breathable privacy cover hides gluteal shelf.

Push to lock brakes integrated into side guard for accessibility.

Ergonomic backrest for maximum handrim access.



INBUILT TRANSPORTATION SOLUTION:

Fold down backrest

Lift legrest section

Tilt wheelchair onto transportation bearings

Easily move wheelchair through doorways or narrow spaces



| Overall Width | Seat Width | Seat Depth | Seat Height | Backrest Height | Total Length | Product Weight | S.W.L. |
|-----------------------------|---------------------------------|------------|---------------------------------|-----------------|--------------|----------------|--------|
| SW+20cm (SP) SW+14cm (T) | 50.5 - 75.5cm (in 5cm steps) | 46 - 51cm | 38 - 45.5cm (in 2.5cm steps) | 45 - 52.5cm | 108cm | 27kg | 300kg |