

BRODA Midline Positioning Wheelchair

The Midline Positioning Wheelchair provides therapeutic pressure relieving comfort and enables better engagement with your environment. It incorporates tilt and recline features and offers easy manoeuvrability to encourage independence. The Midline can be configured to accommodate each individual's needs. The Comfort Tension Seating system provides pressure redistribution and airflow for increased comfort and support. The Midline has various adjustments to ensure an optimised fit, including height adjustable arms, adjustable removable lateral shoulder supports, flip-down foot supports and angle adjustable footplates.

Broda Comfort Tension Seating Powder coated tubular steel frame Height adjustable armrests Adjustable seat tilt and back recline (with Trendelenburg capabilitly) Adjustable wings for upper lateral support Single flip-down, angle-adjustable footrest Elevating, length adjustable legrest Includes headrest, shoulder bolsters, sides, calf, sole, seat and back pads Additional Positioning Padding packages available Tilt and recline provides safe positioning without the use of restraints. Up to 38° tilt and 90° full back recline, including Trendelenburg position. Provides exceptional positioning and support, decreasing postural deviations.



Seat Width	Seat Depth	Seat Height	Tilt	Recline	Warranty	S.W.L.
46 51 56cm	43 46cm	43 - 56cm	38°	90°	2 years	159kg

WA (08) 9456 1661

VIC (03) 9580 1055

SA (08) 7225 8720 www.hospequip.com.au

QLD (07) 3448 9243

NSW 1300 467 737