

V-Drive

The V-Drive system is designed to take over all the effort of pushing from your pushing assistant. Able to handle hills or inclines, it easily navigates all the bumps in the road and ensures comfortable everyday use.

The simple user-friendly control allows for easy adjustment of the speed to ensure a smooth ride and it can be mounted either on the left or right. It also allows you to use the V-drive to propel the wheelchair forwards or backwards.

If you do not wish to use the V-drive, the drive mechanism can simply be lifted up and hung on the push handle until needed again.

To dismantle to fit in the car, it is a simple matter of disconnecting the driving system and battery pack and you're done.



Compatible with a wide variety of wheelchair makes and models

Ideal for longer journeys and hilly terrain

Quick and simple installation and removal

Control can be placed on left or right

Assists in forward and reverse motion

Can be tilted up out of the way if not required

Comes with a convenient transport bag



Controller allows you to adjust speed, select forward or backward propulsion and monitor battery level

	Seat Width	Max speed	Max incline	Product Weight	Battery	Range	S.W.L.
Standard	39 - 48cm	6km/h	6°	19.7kg	2 x 12V / 12Ah	15km	135kg
H.D.	50 - 60cm	6km/h	6°	23kg	2 x 12V / 12Ah	15km	200kg