

# Alber E-Fix



With the e-fix electrical add-on drive for your wheelchair, you can cover longer distances quicker and more easily.

The simple conversion involves swapping the drive wheels and installing the battery pack and control unit. That's it. This means you can continue to use your customised optimally adjusted wheelchair. It is also possible to uncouple the drive wheels to allow for normal manual propulsion.

The quiet, yet powerful drive units are ideal for both indoors and outdoors. Being able to vary your speed from 0.6 to 6km/h allows you more manoeuvrability in confined spaces yet provides enough power to safely negotiate slopes. With a battery that weighs only 2kg, the lightweight slimline design ensures that installation and removal is easy with even one hand and most folding wheelchairs can even be folded up with the battery in place.

Quick release wheels - easy attachment and removal for transport

Configurable - can be individually programmed to suit personal driving characteristics

Durable, lightweight Li-Ion battery with 2 year warranty

EasyConnex magnetic plug system makes connecting the control unit and battery charger even easier

Swivel arm for the controller makes it possible to drive up close to tables



Ergonomic, intuitive controller with large display and optional attendant control makes the e-fix the ideal push and brake aid.



Lightweight components easily attach to your existing wheelchair for a new degree of mobility.

	Wheel Size	Speed	Product Weight	Li-Ion Battery	Range	S.W.L. (Max User)
E35	22"   24"	0.5 - 6km/h	7.9kg (per wheel)	36V 6.0Ah	Up to 16km	170kg (120kg)
E36	24"	0.5 - 6km/h	7.8kg (per wheel)	36V 7.5Ah	Up to 18km	210kg (160kg)