

Alternating Mattresses



BetterLiving alternating pressure surfaces are designed to offer a cost-effective solution for the prevention and management of pressure related injuries.



MULTI-MODE ALARM



QUICK RELEASE CPR



AUTOAIR



2 YEAR WARRANTY

BetterLiving Mattress Replacement / Overlay Systems

Harnessing the latest technical features for improved clinical care and ease of use, BetterLiving offers both a mattress replacement and a mattress overlay option. Ideally suited to home care or aged care environments, these products can be used for people at risk to high risk of pressure injury.

BetterLiving is a brand of cost effective yet clinically effective products designed for the Australian healthcare market. Our extensive product range includes beds and accessories, bathroom and hygiene, seating and positioning, mobility equipment and a range of aids for daily living.

As specialists in the supply of pressure care and manual transfer equipment solutions, Novis is pleased to launch our first range of affordable alternating pressure surfaces by BetterLiving.

Our focus is to offer the most efficient products at the most cost effective price – without compromise to product integrity or clinical patient care.



MATTRESS OVERLAY

MATTRESS REPLACEMENT

SAFETY

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that delivers fail safe patient comfort levels irrespective of their position.

With the AutoCair feature active, the system adjusts to the appropriate pressure setting based on the patient's weight, reducing the need for carer training or monitoring. Manual adjustments can still be made for individualised therapy or comfort.

QUALITY

While we have minimised costs to bring you the most affordable products on the market, we have never minimised on quality.

Our BetterLiving pressure surfaces conform to all the relevant Australian standards. Made using the best quality materials, every product must pass our stringent Quality Control inspection and regular testing protocols.

SIMPLICITY

Offering all the features you need, the BetterLiving pressure surfaces are designed for basic ease of use and minimal maintenance. Easy to set up, simple to operate and straightforward to clean and maintain, these products require minimal training or support.

WATER RESISTANT, MULTISTRETCH COVER

Durable, latex free, multi-stretch vapour permeable cover. Seams are welded and sealed for maximum infection control.

ALTERNATING CELLS

The Mattress Replacement features pure TPU material for maximum immersion and cell-in-cell design for bottoming out protection. The Overlay features lightweight nylon cells.

CPR 2

Quick release CPR access on the Control Unit makes it easy to rapidly deflate the mattress for emergency CPR requirements.



CPR 1

Quick release CPR access on the mattress makes it easy to rapidly deflate the mattress for emergency CPR requirements.

CORD RETENTION

Zippered cord retention sleeve to run mains power cord. Concealed under the waterfall skirt, this feature reduces trip hazards and falls risk.

HANDLES

Sturdy and ergonomic woven handles allow for movement of the inflated mattress. With two handles per side, options for carrying are maximised for the carer or user.

FIXING STRAPS

Easy to use buckle straps secure the mattress to the bed platform for improved profiling and positioning.

DIGITAL CONTROL UNIT

Quick connect umbilical reduces the occurrence of damaged fittings. The upward facing user interface offers more ergonomic access to functions and settings.



Alternating Mattresses

ITEM		Mattress Replacement		Mattress Overlay	
SIZE		Single	King Single	Single	King Single
CODE		APMBL-R02	APMBL-R02K	APMBL-L02	APMBL-L02K
DIMENSIONS	WIDTH	900 mm	1050 mm	900 mm	1050 mm
	LENGTH	2030 mm		2030 mm	
	HEIGHT	200 mm		130 mm	
CAPACITY		200 kg		160 kg	
FEATURES		<ul style="list-style-type: none"> ~ Essential features for pressure injury prevention and treatment ~ Suited for patients at risk to high risk ~ Replace existing mattress for superior support ~ Cell in cell for enhanced user safety and comfort ~ Automated care and therapy ~ Easy to operate and maintain 		<ul style="list-style-type: none"> ~ Feature rich overlay for pressure injury prevention ~ Suited for patients at risk ~ Used over existing mattress for basic support ~ Automated care and therapy ~ Easy to operate and maintain 	
CONFIRMS		ARTG 289458, IEC 60601-1, IEC60601-1-2, IEC60601-1-11			
WARRANTY		2 year			



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.