



 **Bodypoint**[®]

STRENGTHEN YOUR POSITION



Josh Dueck Professional Skier, Paralympian and Bodypoint® User
Photo: Alain Sleighter

4-11

PELVIC POSITIONING

- Hip Belts
- Sub-ASIS Pads and Compatible Belts
- Evoflex® Pelvic Stabilizers
- Hip Belt Accessories
- Slip-On Pads
- Leg Harnesses

12-15

LOWER BODY POSITIONING

- Ankle Huggers® Support Straps
- Adjustable Angle Footplates
- Calf Panels and Calf Straps

16-22

UPPER BODY POSITIONING

- Stayflex™ Chest Supports
- Shoulder Harnesses
- Aeromesh® Rapid-Dry Bath Belts
- Monoflex™ Chest Supports
- Chest Belts
- PivotFit™ Shoulder Harnesses

23-25

POWERCHAIR COMPONENTS

- Midline Joystick Mounting System
- Joystick Handles

26-27

HARDWARE

- Fixit™ and Clickit™
- FrameSaver™ Clamps
- Frame Mounting Clamps
- Quick-Release Connectors
- Seat Tube Mounting Bracket Kit
- Belt Mounting Extension Kit

28

WARRANTY AND TRADEMARK INFORMATION



Standards-tested products can be trusted to deliver the fit, durability, safety and performance you expect and deserve. As a result, you can have greater confidence in the devices you sell, buy, apply and wear—empowering you to be a positioning champion.



Standards help strengthen your position

“There’s just no one-size-fits-all solution for achieving individual potential. At Bodypoint® we understand and embrace this challenge.”

—Dave Hintzman
President and Co-Founder,
Bodypoint®

Since our inception, Bodypoint has been the leading voice in the development of global standards, driving improvements in the design, quality and application of postural support devices.

After hard work, our efforts paid off!

The **International Standards Organization (ISO)**, the **American National Standards Institute (ANSI)** and the **Rehabilitation Engineering and Assistive Technology Society of North America (RESNA)** have introduced new standards for wheelchair seating.

These new standards not only provide sizing

guidelines and performance criteria, they also reflect a growing recognition across the wheelchair community that raising the bar—and our expectations—can empower us all to achieve our full potential.



Meet **Barry Long**, featured above, or other satisfied customers at:

www.bodypoint.com/success-stories.aspx



Learn more today at
www.bodypoint.com/standards.aspx



PELVIC POSITIONING

The pelvis is the key to postural alignment as it dictates the position of the trunk, head and extremities.



PUSH-BUTTON BUCKLE

The Bodypoint® Push-Button Buckle is sleek and low profile. Engineered for comfort and performance, this new design offers a range of features that therapists and users have requested for years. The patented construction allows for removal of the buckle cover for easy cleaning, or to change the desired button hole size – in effect providing three buckle options in a single design.

Light Weight • Durable • Low Profile • Easy to Clean • Assisted Spring Release



Lighter weight (1/3 lighter than comparable push-button models) with rounded edges that improve comfort and prevent snagging on clothing.



Chrome-plated steel construction for superior strength and durability. Lower profile for a less bulky appearance.



Optional covers:
Reduced Access Cover (Ø14mm button),
Security Cover (Ø6mm button).
Sold separately.
See page 11 for ordering information.

SELECT A BODYPOINT HIP BELT IN 6 STEPS

1. SELECT BELT TYPE

- Padded (pg. 7)
- Non-Padded (pg. 11)
- Hook-and-Loop Compatible (pg. 8)

2. TWO- OR FOUR-POINT



TWO-POINT HIP BELT

Two points of attachment to the wheelchair. Suitable when the belt is supporting users with low tone or weakness.



FOUR-POINT HIP BELT

Four points of attachment to the wheelchair. Primary straps position the belt like two-point hip belt and secondary straps anchor the belt in position to prevent it from riding up into the abdomen or twisting.

3. SELECT BUCKLE

OPTIONAL COVERS:



Push-Button Buckle

Low profile, smooth, snag-free corners. Removable cover for changing release hole size and easy cleaning. Standard Access Cover (a) included. Reduced Access Cover (b) and Security Cover (c) sold separately.



Plastic Side-Release
Lightweight, secure snap action.



Rehab Latch™ Buckle
Lightweight, strong, low profile, easy to operate with limited hand function.



Extra Small Metal Push-Button
Strong, secure, child-sized buckle. Buckle Security Cover also available. (See pg. 11)

4. SELECT A PULL TYPE



CENTER-PULL HIP BELT

Operator: User or caregiver
Adjustment: Tightens at the center buckle.
Buckles: Available in the full buckle range.
Applications: Clients with low muscle tone or weakness.



REAR-PULL HIP BELT

Operator: Caregiver
Adjustment: Tightens from the rear of the pad pulling toward the center.
Buckles: Available in full buckle range.
Applications: Clients with high tone. Counters pelvic rotation by allowing the caregiver to position one side of the pelvis, lock it in place and then position on the other side.



DUAL-PULL HIP BELT

Operator: User or caregiver
Adjustment: Tightens from both sides of the center buckle.
Buckles: Side-release buckle only
Applications: Clients requiring greater adjustment to accommodate clothing and weight changes.

5. SELECT A SIZE

Sizes range from extra small to large and vary by style (see pg. 7-8). See measurement section to the left for details.

6. SELECT A MOUNTING OPTION




Cinch-Mount™
Time saving, precise, easy adjustment

Flat-Mount™
Versatile, lowest profile




FrameSaver™ Clamp (optional)
Create a mounting point without drilling. Fits Ø7/8"–Ø1" (22-25mm) tube.

MEASURING FOR A HIP BELT



Measure hip width across the greater trochanters with the person seated. If the hip width falls between sizes on our chart, consider other factors such as growth, weight changes and clothing.

XS	5–9" (13–23cm)
S	7–11" (18–28cm)
M	9–15" (23–38cm)
L	13–19" (33–48cm)



	A	B	C	D
XS	1" (25mm)	1-3/4" (44mm)	6" (15cm)	50" (127cm)
S	1-1/2" (38mm)	2-1/4" (57mm)	7" (18cm)	55" (140cm)
M	1-1/2" (38mm)	2-1/2" (64mm)	9" (23cm)	60" (152cm)
L	2" (50mm)	3" (76mm)	11" (28cm)	65" (165cm)



All Bodypoint® hip belts are tested to the new standards for wheelchair seating postural support devices.

Standards-tested products can be trusted to deliver the fit, durability, safety and performance you expect and deserve—empowering wheelchair users, care providers and commercial partners to be positioning champions.

TWO-POINT PADDED HIP BELTS

		Push-Button	Plastic Side-Release	Rehab Latch™	Extra Small Metal Push-Button
CENTER-PULL	XS		HB215XS		HB235XS†
	S	HB205S*	HB215S	HB225S	
	M	HB205M*	HB215M	HB225M	
	L	HB205L**	HB215L	HB225L	
DUAL-PULL	XS		HB216XS		
	S		HB216S		
	M		HB216M		
	L		HB216L		
REAR-PULL	XS		HB219XS		HB239XS†
	S	HB209S*	HB219S	HB229S	
	M	HB209M*	HB219M	HB229M	

*Accepts the FS018D-1 Reduced Access and FS018C-1 Security Covers (see page 11).

**Accepts the FS019D-1 Reduced Access and FS019C-1 Security Covers (see page 11).

†Accepts FS029C-1 Buckle Security Cover (see page 11).

FOUR-POINT PADDED HIP BELTS

		Push-Button	Plastic Side-Release	Rehab Latch™	Extra Small Metal Push-Button
CENTER-PULL	XS		HB415XS		HB435XS†
	S	HB405S*	HB415S	HB425S	
	M	HB405M*	HB415M	HB425M	
	L	HB405L**	HB415L	HB425L	
DUAL-PULL	XS		HB416XS		
	S		HB416S		
	M		HB416M		
	L		HB416L		
REAR-PULL	XS		HB419XS		HB439XS†
	S	HB409S*	HB419S	HB429S	
	M	HB409M*	HB419M	HB429M	

*Accepts the FS018D-1 Reduced Access and FS018C-1 Security Covers (see page 11).

**Accepts the FS019D-1 Reduced Access and FS019C-1 Security Covers (see page 11).

†Accepts FS029C-1 Buckle Security Cover (see page 11).


To complete your part number, select the mounting options below:

Flat-Mount™ + FrameSaver™ Clamp		-A1
Flat-Mount™ Only		-B1
Cinch-Mount™ + FrameSaver™ Clamp		-A2
Cinch-Mount™ Only		-B2
FrameSaver™ Clamp, Cinch-Mount™ and Flat-Mount™ hardware is standards-tested as part of the complete hip belt system.		

Example: **HB205S-B1**

Product number for the Center-Pull Push-Button Buckle, size Small with Flat-Mount™ only.

SUB-ASIS COMPATIBLE BELTS



SUB-ASIS PADS

The teardrop-shaped gel pads comfortably space the belt away from the abdomen (reducing pressure on the bladder), while firmly controlling the pelvis just below the ASIS.

PD402-2 Extra Small/Small

PD401-2 Medium/Large



To complete your part number, select the end-fitting option below:

Flat-Mount™ + FrameSaver™ Clamp		-A1
Flat-Mount™ Only		-B1
Cinch-Mount™ + FrameSaver™ Clamp		-A2
Cinch-Mount™ Only		-B2

Combined with Sub-ASIS pads, our Sub-ASIS Compatible Belt provides two firm points of contact with the pelvis. The result is greatly enhanced stability and improved control of obliquity and rotation while minimizing abdominal compression.

		Size	Push-Button	Plastic Side-Release	Rehab Latch™	Extra Small Metal Push-Button	
TWO-POINT	Center-Pull	XS				LB235XS ¹	
		M	LB205M [*]		LB225M		
		L	LB205L ^{**}		LB225L		
	Rear-Pull	M			LB229M		
FOUR-POINT	Center-Pull	XS				LB435XS ¹	
		S			LB425S		
		M	LB405M [*]		LB425M		
		Dual-Pull	L	LB405L ^{**}		LB425L	
		Rear-Pull	XS		LB416XS		
		M			LB429M		

^{*}Accepts the FS018D-1 Reduced Access and FS018C-1 Security Covers (see page 11). ^{**}Accepts the FS019D-1 Reduced Access and FS019C-1 Security Covers (see page 11). ¹Accepts FS029C-1 Buckle Security Cover (see page 11).

Measuring for Sub-ASIS Compatible Belts

Size	Sitting Hip Width*	User Mass
XS	8-11" (20-28 cm)	15-35 kg
M	11-14" (28-36 cm)	35-70 kg
L	14-18" (36-46 cm)	70-130 kg

*Sitting hip width measured across the greater trochanters



All Bodypoint® hip belts are tested to the new standards for wheelchair seating postural support devices.

Standards-tested products can be trusted to deliver the fit, durability, safety and performance you expect and deserve—empowering wheelchair users, care providers and commercial partners to be positioning champions.

THE EVOFLEX® PELVIC STABILIZER

Takes the Hassles out of Your Hip Belt



The Evoflex® Pelvic Stabilizer is the first support with stiffened end straps that stay right where you need them and won't twist or fall into the wheels.

- Better stability due to a closer fit
- Resistant to abrasion, tearing and temperature extremes

Other features include:

Slim, tapered strap design improves stability by allowing a closer fit

Unique, stiffened end straps keep belt from twisting or falling into the wheels

Simple, low-profile fastener fits into even the tightest places

See Hardware section for more information on these mounting options.



Pivot belt either forward or back for easier transfers.




Mounting to manual wheelchair frame using Band Clamps.



Mounting to the seat rail using Seat Tube Mounting Bracket Kit.

MEASURING FOR AN EVOFLEX® PELVIC STABILIZER

- 1** Measure seated hip width of the user across the greater trochanters, then determine your mounting points on the chair.



RECOMMENDED FIT RANGE

Size	Sitting Hip Width	Weight
XS	20-28cm (8-11")	15-35kg (33-77lbs)
M	28-36cm (11-14")	35-70kg (77-154lbs)
L	36-46cm (14-18")	70-130kg (154-298lbs)

EVOFLEX® PELVIC STABILIZER

PUSH-BUTTON

Center-Pull, Two-Point

EB205L	Large
EB205M	Medium

REHAB LATCH™


Center-Pull, Two-Point

EB225L	Large
EB225M	Medium

METAL PUSH-BUTTON

Center-Pull, Two-Point

EB235XS	X-Small
---------	---------



- 2** Select the preferred mounting point on the chair or seating system, and measure between them across the user to determine the needed belt length. You may want to consider other factors such as seat cushion thickness*, growth, and weight changes/clothing.

Size	A Strap Width	B Maximum Mounting Length	C Minimum Mounting Length
XS	2.5cm (1")	63cm (24-3/4")	43cm (16-3/4")
M	4cm (1-1/2")	76cm (29-3/4")	50cm (19-3/4")
L	5cm (2")	96cm (37-3/4")	65cm (25-3/4")

* See Hardware section for Belt Mounting Extension Kit option.



BARRY LONG

Woodinville, WA, USA

Evoflex® User Profile

Barry's unbridled enthusiasm for life and adventure is contagious. An intrepid explorer, his extreme sports exploits give him a unique way to test the limits of his equipment. We asked Barry, who is not one to mince words, to product test the Evoflex® Pelvic Stabilizer and provide us with feedback.



See what Barry has to say at

[www.bodypoint.com/
success-stories-barry-long.aspx](http://www.bodypoint.com/success-stories-barry-long.aspx)

**OPTIONAL PUSH-BUTTON
BUCKLE COVERS**

(Push-Button Buckle shown on page 6)



REDUCED ACCESS COVER

FS018D-1 1-1/2", Ø14mm button

For use with small and medium belts.

FS018D-50 Pack of 50

FS019D-1 2", Ø14mm button

For use with large belts.

FS019D-50 Pack of 50



SECURITY COVER

FS018C-1 1-1/2", Ø6mm button

For use with small and medium belts.

FS018C-50 Pack of 50

FS019C-1 2", Ø6mm button


For use with large belts.

FS019C-50 Pack of 50



**EXTRA SMALL
BUCKLE SECURITY COVER**

Snaps securely onto our Extra Small Metal Push-Button Buckle, limiting access by even the smallest fingers. Open instantly with a pen tip. **FS029C-1**



Shown: Previous Metal Push-Button Buckle with Buckle Security Cover Attached.



BUCKLE SECURITY COVER

Snaps securely onto our previous Metal Push-Button Buckle (no longer available), limiting access by even the smallest fingers. Open instantly with a pen tip. **FS025C-1**

**NON-PADDED
HIP BELTS**



	Size	Push-Button	Plastic Side-Release	Rehab Latch™	Compact Push-Button
Center-Pull	XS		NB215XS		
	M	NB205M*	NB215M	NB225M	NB245M
	L	NB205L**	NB215L	NB225L	NB245L'
Center-Pull Quad Belt	L			NB225L-T	

*Accepts the FS018D-1 Reduced Access and FS018C-1 Security Covers.
**Accepts the FS019D-1 Reduced Access and FS019C-1 Security Covers.

To complete your part number, select the end-fitting option in the chart below



Flat-Mount™ + FrameSaver™ Clamp		-A1
Flat-Mount™ Only		-B1
Cinch-Mount™ + FrameSaver™ Clamp		-A2
Cinch-Mount™ Only		-B2

SLIP-ON PADS

Allow you to add cushioning where additional relief is needed on belts and harnesses.

PD304 Small

Fits webbing: 1" (25mm)
Length: 8" (20cm)

PD305 Medium

Fits webbing: 1-1/2" (38mm)
Length: 10" (25cm)

PD306 Large

Fits webbing: 2" (50mm)
Length: 12" (30cm)



LEG HARNESS

Holds each leg at the inner thigh, providing exceptional resistance to pelvic rotation and obliquity, without limiting anterior movement or pressing on the bladder. Sold in pairs.


LH600S-A1 or -A2 or -B2 Small
Length: 9" (23cm), Width: 1-1/2" (4cm)

LH600M-A1 or -A2 or -B2 Medium
Length: 13" (33cm), Width: 1-3/4" (4.5cm)

LH600L-A1 or -A2 or -B2 Large
Length: 15" (38cm), Width: 2" (5cm)



**MEASURING FOR
A LEG HARNESS**



Measure from the Anterior Superior Iliac Spine (ASIS) to the point where the inner thigh touches the seat.

S Fits 8-11" (20-28cm)

M Fits 11-14" (28-36cm)

L Fits 14-16" (36-41cm)

LOWER BODY POSITIONING

Lower body positioning affects the position of the pelvis, balance and freedom of movement of the upper body.





Hook-and-Loop Closure

ANKLE HUGGERS

SUPPORT STRAPS

Ankle Huggers® Support Straps are the original patented solution for stabilizing the feet while allowing controlled movement. They reduce joint stress by holding above the ankle, and greatly improve comfort and tolerance for positioning. Engineered for outstanding fit and lasting quality with our patented Laminar Pad™ and proprietary webbing. Ankle Huggers Support Straps are the most secure means available for holding the foot.



Patented reinforcement patch provides enhanced strength between the straps and pads.



Circumferential pad shape provides even pressure distribution.



Side-release attachment for exact fit. Also available in hook-and-loop closure for easy adjustment.



**HOOK-AND-LOOP
RELEASE ANKLE HUGGERS®
SUPPORT STRAPS**

Sold in pairs (1 right/1 left)
Soft hook-and-loop adjustment strap conforms to ankle curvature and is quickly tightened with each use to ensure a snug fit.


- FT210XS** Extra Small
- FT210S** Small
- FT210M** Medium
- FT210L** Large
- FT210XL** Extra Large

**SIDE-RELEASE ANKLE HUGGERS®
SUPPORT STRAPS**

Sold in pairs (1 right/1 left)
Adjustment strap retains its precise position with every use, and side-release buckle ensures strong closure even with extensive use.

- FT220XS** Extra Small
- FT220S** Small
- FT220M** Medium
- FT220L** Large
- FT220XL** Extra Large

**MEASURING
FOR ANKLE HUGGERS®
SUPPORT STRAPS**



Measure circumference above ankle.

XS	5-1/2"-7" (14cm-17cm)
S	6-1/2"-8" (17cm-20cm)
M	7-1/2"-9" (19cm-23cm)
L	8-1/2"-11-1/2" (22cm-29cm)
XL	11"-13" (28cm-33cm)



AEROMESH® CALF PANEL

Features our proprietary Aeromesh® material. Aeromesh® fabric is bonded to a durable stiffener to allow air flow while maintaining shape.

- SP103S** Small
- SP103M** Medium
- SP103L** Large



AEROMESH® PADDED CALF STRAPS

Calf Straps provide comfortable support, keeping feet safely on the footplates. Retains its shape for the best pressure distribution and durability.

- SP102S** Small
- SP102M** Medium
- SP102L** Large




HOOK-AND-LOOP TOE STRAPS

Toe strap slots are provided on all Bodypoint footplates to allow for easy attachment. Sold in pairs.


FT221

MEASURING FOR CALF PANELS AND CALF STRAPS




For proper fit, measure chair width.

S	12"-15" (31cm-38cm)
M	15"-19" (38cm-48cm)
L	19"-22" (48cm-56cm)



Calf Strap Dimensions:

	A	B	C
S	31-1/2" (80cm)	14" (36cm)	3-1/2" (89mm)
M	35-1/2" (90cm)	16" (41cm)	3-1/2" (89mm)
L	39-1/2" (100cm)	18" (46cm)	3-1/2" (89mm)



Calf Panel Dimensions:

	A	B	C
S	31-1/2" (80cm)	14" (36cm)	7-1/2" (19cm)
M	35-1/2" (90cm)	16" (41cm)	8-1/4" (21cm)
L	39-1/2" (100cm)	18" (46cm)	9-1/4" (23cm)

FOOTPLATE MOUNTING OPTIONS



Mounted in hanger tubes



Mounted high for short leg lengths (requires HW404-2 or HW406-2 clamp kits)



Mounted at different heights for single leg amputee (requires HW404-2 or HW406-2 clamp kits)



L-SHAPED EXTENSION TUBE
Ø3/4" (19mm)

FT021 Length: 8-1/2" (22cm)
FT022 Length: 12-1/2" (32cm)



FULCRUM SERIES EXTENSION TUBE
Ø3/4" (19mm)

FT023 Length: 8-1/2" (22cm)
FT024 Length: 12-1/2" (32cm)



MOUNTING CLAMP
HW404 Ø3/4" (19mm) X Ø7/8" (22mm)
HW404-2 Pair
HW406 Ø3/4" (19mm) X 1" (25mm)
HW406-2 Pair



FOOTPLATE CLAMP
Clamp permits easy adjustment of footplate and locks tightly in position for absolute stability.
FT032 Ø3/4" (19mm)



FOOTPLATE
Heavy gauge aluminum with super-tough, non-skid nylon surface stands up to hard use.
FT051 Small
Width: 4" (10cm)
Length: 7-3/4" (20cm)
FT051-2 Pair
FT052 Large
Width: 4-3/4" (12cm)
Length: 9-1/4" (24cm)
FT052-2 Pair

ADJUSTABLE-ANGLE FOOTPLATE KITS

Includes: choice of extension tubes (pair), choice of footplates (pair), and footplate clamps (pair).

Footplate Size	L-SHAPED		Footplate Size	FULCRUM, Ø3/4" (19mm) TUBE	
	Tube Length			Tube Length	
	8-1/2" (22cm)	12-1/2" (32cm)		8-1/2" (22cm)	12-1/2" (32cm)
S	FT101	FT101L	S	FT105	FT105L
L	FT102	FT102L	L	FT106	FT106



ADJUSTABLE ANGLE FOOTPLATE
FT111 Small
FT111-2 Pair
FT112 Large
FT112-2 Pair



UPPER BODY POSITIONING

Trunk support is essential for good respiration, head positioning and functional reach with stability.



All Bodypoint® upper body positioning products are tested to the new standards for wheelchair seating postural support devices.

Standards-tested products can be trusted to deliver the fit, durability, safety and performance you expect and deserve—empowering wheelchair users, care providers and commercial partners to be positioning champions.



STAYFLEX CHEST SUPPORT

The patented design of the Bodypoint® Stayflex™ includes a stabilizing lower panel to limit stretching and riding up. The proprietary elastic pad material allows a comfortable, controlled amount of stretch. Plus, the lower quick-release swivel buckles self-adjust to the line of tension for proper fit and placement.



Proprietary elastic pad material offers controlled stretch for dynamic movement.



Rear-pull attachment is ideal for use with thicker backrests or when you want the buckle or adjustment away from the user.




Stabilizing lower panel limits stretching and riding up up to reduce choking hazard.

STANDARD STAYFLEX™		FRONT-PULL		REAR-PULL	
 Shown with zipper	Size	With Zipper	Without Zipper	With Zipper	Without Zipper
	XS		SH311XS		SH321XS
	S	SH310S	SH311S	SH320S	SH321S
	M	SH310M	SH311M	SH320M	SH321M
	L	SH310L	SH311L	SH320L	SH321L
	XL	SH310XL	SH311XL	SH320XL	SH321XL

NARROW STAYFLEX™		FRONT-PULL		REAR-PULL	
 Shown with zipper	Size	With Zipper	Without Zipper	With Zipper	Without Zipper
	M	SH330M	SH331M	SH340M	SH341M
	L	SH330L	SH331L	SH340L	SH341L
	XL	SH330XL	SH331XL	SH340XL	SH341XL

MEASURING FOR STAYFLEX™ AND SHOULDER HARNESS




Measure shoulder width. If the shoulder width falls between two sizes, consider other factors such as growth, weight changes, and clothing.

XS	9-1/2-11" (24-28cm)
S	11-13" (28-33cm)
M	13-16" (33-41cm)
L	16-19" (41-48cm)
XL	19-22" (48-56cm)


Choose the mounting option best suited for your application. Specify the Flat-Mount™ option (-B1) when looking for the best value or the Cinch-Mount™ option (-B2) when ease of installation is important.

STAYFLEX™ MEASUREMENTS



STANDARD STAYFLEX™


Size	A	B	C	D	E
XS	8-1/2" (22cm)	5" (13cm)	5-1/4" (13cm)	1-1/2" (4cm)	3-3/4" (9cm)
S	9" (23cm)	5-1/2" (14cm)	6-1/2" (16cm)	1-3/4" (4cm)	4-3/4" (12cm)
M	12" (30cm)	7" (18cm)	8" (20cm)	1-3/4" (4cm)	5" (13cm)
L	15" (38cm)	8-1/2" (21cm)	9" (23cm)	2-1/4" (6cm)	5-1/4" (13cm)
XL	18" (46cm)	10" (25cm)	10-1/2" (27cm)	2-1/4" (6cm)	6-1/2" (16cm)



NARROW STAYFLEX™


Size	A	B	C	D	E
M	12" (30.5cm)	7" (18cm)	8" (20.5cm)	1-3/4" (4cm)	3-1/2" (9cm)
L	15" (38cm)	8-1/2" (21.5cm)	9" (23cm)	2-1/4" (6cm)	3-1/2" (10cm)
XL	18" (46cm)	10" (25.5cm)	10-1/2" (27cm)	2-1/4" (6cm)	3-3/4" (10cm)

STAYFLEX™ AND SHOULDER HARNESS ATTACHMENT OPTIONS



REAR-PULL

Rear-pull attachment is ideal for use with thicker backrests or when you want the buckle or adjustment away from the user.



OPTIONAL CONFIGURATION

The included side-release buckles may be removed from the rear-pull strap to provide a continuous length of straight webbing for use with cam buckles or other mounting hardware.

Cam buckles shown not included. See page 26 for ordering information.



FRONT-PULL

Front-pull attachment provides the maximum stability while allowing adjustment standing either behind or in front of the user.

TRIMLINE-STYLE SHOULDER HARNESS Made of two independent padded straps which allow the greatest degree of freedom in fitting the chest.

	Size	Pad Length	Front-Pull	Rear-Pull
	S	9-1/4" (23.5cm)	SH230S	SH240S
	M	11-1/2" (29cm)	SH230M	SH240M
	L	13-1/4" (34cm)	SH230L	SH240L



CROSSOVER POSITION
Comfortable fit which does not restrict arm movement.

BACKPACK POSITION

Promotes respiration by encouraging shoulder retraction while stabilizing the trunk.

AEROMESH® RAPID-DRY BATH BELT



FEEL SAFE AND SECURE with the Aeromesh® Rapid-Dry Bath Belt



AEROMESH® RAPID-DRY BATH BELT

The Bodypoint bath belt with our proprietary Aeromesh® padding provides safety without compromising the comfort and health of the user. Water and air easily pass through the Aeromesh fabric, allowing the client to be bathed completely and the belt to dry quickly after use. Its flexible construction and universal shape allow it to be used in various ways to secure the user: across the chest or abdomen to prevent forward/lateral falls, or fitted snugly across the knees to keep clients from slipping forward and out of the chair.

Two-piece Bath Belt

BB101M-1 Medium (without clamps)

BB101L-1 Large (without clamps)

One-piece Bath Belt


BB111L-1 Large, grommets

BB111XL-1 X-Large, grommets

BB112L-1 Large, plastic slides

BB112XL-1 X-Large, plastic slides

Water and air easily pass through the Aeromesh® fabric so it dries quickly.



XL UNIVERSAL ELASTIC STRAP AND ELASTIC EXTENSION STRAPS



Elastic Extension Straps

The Elastic Extension Straps were created to give added length to the original Universal Elastic Strap while keeping all the same great features. This accessory simply attaches to the end of the Universal Elastic Strap using the hook- and-loop closure to extend the overall length. The thumb loop end attaches to the Elastic Extension Strap to secure the end. Usable Length: 72.5cm (28-1/2")

Universal Elastic Straps

For temporary positioning in the therapy clinic, during fittings, when rolling over rough streets, or for sports, the Universal Elastic Strap is the stable solution. This strong elastic strap can be used on the chest, abdomen, lap or legs, and has an easy-grasp end tab with a thumbhole for users with limited hand function – truly a versatile support and positioning tool for caregivers, therapists and users. Usable Length: 142cm (56")

- Durable, high-stretch neoprene with soft nylon loop fabric
- Simple hook-and-loop closure (similar to Velcro®) for easy one-handed opening and adjustment
- Easily cut-to-length without fraying



The Universal Elastic Straps are sold individually and are now in X-Large:

SP110S Small, width: 8cm (3")
 SP110L Large, width: 11cm (4-1/2")
 SP110XL X-Large, width: 15cm (6")

The Elastic Extension Straps are sold individually and are available in all corresponding sizes for the Universal Elastic Straps:

SP111S Small, width: 8cm (3")
 SP111L Large, width: 11cm (4-1/2")
 SP111XL X-Large, width: 15cm (6")



MONOFLEX™ CHEST SUPPORT

The perfect dynamic abdominal or chest belt. Comfortable, high-strength elastic pad provides just the right amount of stretch while maintaining substantial control. Side-release buckles allow adjustment and release on both sides for convenience. Versatile mounting hardware fits any type of solid or sling backrest.

Size	Pad Length	Pad Width	Part No.
XS	12-1/4" (31cm)	3" (8cm)	SH103XS
S	15-1/2" (39cm)	3-1/2" (9cm)	SH103S
M	18-1/4" (46cm)	4" (10cm)	SH103M
L	21-1/2" (55cm)	4-1/2" (12cm)	SH103L
XL	24-1/2" (62cm)	5" (13cm)	SH103XL



CHEST BELT

Features a premium hook-and-loop closure with a leather-stiffened ring for easy adjustment by users with limited hand function.


Usable Length: 46" (117cm)
 Webbing Width: 2" (50mm)

SH102-B2 Cinch-Mount™
SH102P-B2 Cinch-Mount™
 (with PD306 Slip-On Pad. See pg. 11)

SH102-B1 Flat-Mount™
SH102P-B1 Flat-Mount™
 (with PD306 Slip-On Pad. See pg. 11)

MEASURING FOR MONOFLEX™

Measure the user's trunk around the front and sides at the height at which the Monoflex will be worn. Consult the pad sizes (above) and select the size that will cover the front and wrap at least halfway around each side of the trunk, leaving a gap in the front of the wheelchair's backrest to allow space for tightening the buckles.



Size	Min/Max Fit Range
XS	14-17-1/2" (36-45cm)
S	18-21-1/2" (46-55cm)
M	22-25-1/2" (56-65cm)
L	26-29-1/2" (66-75cm)
XL	30-33-1/2" (76-85cm)

Changing the Shape of Support

PIVOTFIT™ SHOULDER HARNESS



“The PivotFit™ helps me sit up straighter because it holds my shoulders back, so I’m not all flopped over.”

– Jordan F.
Duvall, WA, USA



More on the PivotFit™ and how it’s helped Jordan at www.bodypoint.com/pivotfit-shoulder-harness.aspx

PIVOTFIT™ SHOULDER HARNESS

The form-fitting shape of the PivotFit™ Shoulder Harness offers optimal support and a comfortable, individualized fit. Curving Laminar Pads™ follow the contours of the shoulders and ribcage with a lowered sternum adjustment strap to reduce the risk of choking. Complete with a new buckle and refined sewing construction for easy-care, the PivotFit™ presents a sleek, form-fitting look. Notable features include:

- Easy-cinching sternum strap for quick adjustments over bulky sweaters and jackets
- 25% greater pressure distribution than a traditional H-Style Harness

AVAILABLE IN TWO STYLES:



The **Standard** non-stretch model, with its durable and easy-to-clean polyester covering, reduces movement to address serious positioning needs, like severe kyphosis or trunk rotation.



The **Dynamic** model uses our patented dual-zone technology to enhance comfort by allowing moderate shoulder movement without reducing support at the chest – ideal for helping manage spastic tone.

Patented dual-zone elasticity distributes stretch and pressure evenly for increased comfort and support (Dynamic PivotFit™ only).

Unique shape contours to the shoulders and ribcage for a more supportive, comfortable fit.

New, one-of-a-kind, high-strength swivel buckle allows the pads to pivot into the best placement for the user's body shape and comfort.



Standard PivotFit™ Shoulder Harness

Front-Pull		Rear-Pull	
SH250XL	X-Large	SH260XL	X-Large
SH250L	Large	SH260L	Large
SH250M	Medium	SH260M	Medium
SH250S	Small	SH260S	Small
SH250XS	X-Small	SH260XS	X-Small
SH250XXS	XX-Small	SH260XXS	XX-Small

Dynamic PivotFit™ Shoulder Harness

Front-Pull		Rear-Pull	
SH270XL	X-Large	SH280XL	X-Large
SH270L	Large	SH280L	Large
SH270M	Medium	SH280M	Medium
SH270S	Small	SH280S	Small
SH270XS	X-Small	SH280XS	X-Small
SH270XXS	XX-Small	SH280XXS	XX-Small

Mounting options: Complete your product sku by choosing -B1 Flat-Mount™ or -B2 Cinch-Mount™ (ex. SH250L-B1)



POWER CHAIR COMPONENTS

Improving the accessibility of power wheelchair controls lets the user achieve a higher level of independence.


MIDLINE JOYSTICK MOUNTING SYSTEM

The Midline Joystick Mounting System improves accessibility and function. Control units and communication devices are easily mounted within reach without additional chair modifications. The precision-machined construction withstands the toughest daily use.



MIDLINE ARMS

Midline Arms support joystick boxes and communication devices weighing up to 30lbs.(13.6kg).



PC031 J-Shaped Midline Arm Bracket





PC032 L-Shaped Midline Arm Bracket




PC033 Straight Midline Arm Bracket

MIDLINE JOYSTICK MOUNTING KITS Kits contain: 1 Rotating Shaft, 1 Mounting Clamp, choice of Midline Arm and choice of Mounting Adapter.


KIT STYLE	JOYSTICK MOUNTING ADAPTERS						
	UNIVERSAL	MK6I MPJ+, PSF+, AND PSR+	QUICKIE QTRONIX	MK5 SPJ AND DPJ	PGDT PILOT+	MK5 MPJ, Q-LOGIC NE, NE+, REMOTE, AND EX	PERMOBIL COMPACT JOYSTICKS
	 PC041	 PC042	 PC043	 PC044	 PC045	 PC046	 PC047
J-SHAPED ARM 	PC201	PC202	PC203	PC204	PC205	PC206	PC207
L-SHAPED ARM 	PC201A	PC202A	PC203A	PC204A	PC205A	PC206A	PC207A
STRAIGHT ARM 	PC201B	PC202B	PC203B	PC204B	PC205B	PC206B	PC207B

TRI-LOCK™ ROTATING SHAFT POSITIONS


Shown with a
Straight Midline Arm Bracket.




0° Horizontal driving position




90° Allows user to transfer



90° Press to disengage



240° Moved out of the way for
table access



JOYSTICK HANDLE COLLECTION

Evaluate and select the best handle style for optimal driving control from our collection of nine flex-shaft joystick handles. Choose a display with handles to fit either 3/16" or 1/4" diameter stems.

PC130-A KIT

Fits stem size $\text{\O}3/16''$

PC130-B KIT

Fits stem size $\text{\O}1/4''$

BODYPOINT® HAS JOYSTICK HANDLES to fit every type of user.



3" U-SHAPED

PC101A Fits stem size: $\text{\O}3/16''$

PC101B Fits stem size: $\text{\O}1/4''$



4" U-SHAPED

PC102A Fits stem size: $\text{\O}3/16''$


PC102B Fits stem size: $\text{\O}1/4''$



BALL KNOB

PC103A Fits stem size: $\text{\O}3/16''$

PC103B Fits stem size: $\text{\O}1/4''$



STRAIGHT

PC104A Fits stem size: $\text{\O}3/16''$

PC104B Fits stem size: $\text{\O}1/4''$



RUBBER DOME

PC107A Fits stem size: $\text{\O}3/16''$

PC107B Fits stem size: $\text{\O}1/4''$



2" RUBBER MUSHROOM HEAD

PC108AS Fits stem size: $\text{\O}3/16''$

PC108BS Fits stem size: $\text{\O}1/4''$



2 1/2" RUBBER MUSHROOM HEAD

PC108AL Fits stem size: $\text{\O}3/16''$

PC108BL Fits stem size: $\text{\O}1/4''$



ALUMINUM TOP

PC109A-BLUE Fits stem size: $\text{\O}3/16''$

PC109B-BLUE Fits stem size: $\text{\O}1/4''$

Grey/Red available as -GREY or -RED

HARDWARE

Bodypoint® hardware provides a wide range of mounting options for all your seating equipment.



Bodypoint FrameSaver™ Clamp and Cinch-Mount™ hardware is standards-tested as part of the complete postural support device mounting system.



FRAMESAVER™ CLAMP

Fits Tube Size: Ø7/8" (22mm) and Ø1" (25mm), Optional shim for Ø3/4" (19mm) tube also available (part no. HW303-19SHIM-2 or -100)

- HW303-B-2** Quantity: 2
- HW303-B-10** Quantity: 10
- HW303-B-50** Quantity: 50
- HW303-B-100** Quantity: 100



ALUMINUM FRAME CLAMP

This double screw clamp allows you to separately secure the mounting position prior to attaching your accessory—even on an upholstered seat rail!

- HW310-22-A-2** Ø7/8" (22mm) Quantity: 2
- HW310-22-A-4** Ø7/8" (22mm) Quantity: 4
- HW310-22-A-50** Ø7/8" (22mm) Quantity: 50
- HW310-25-A-2** Ø1" (25mm) Quantity: 2
- HW310-25-A-4** Ø1" (25mm) Quantity: 4
- HW310-25-A-50** Ø1" (25mm) Quantity: 50



QUICK-RELEASE CONNECTOR

The patented Bodypoint Quick-Release Connector provides an easily detached, solid mount that can be quickly released—even with weak hand function. The simple and reliable push-to-close, pull-to-release mechanism can be operated with one hand and locked for extra security when needed. With a slim, 1/2" profile and durable stainless steel and Delrin construction, these versatile mounts fit securely in all but the most inaccessible places.

- HW311** Quantity: 1
- HW311-2** Quantity: 2



CAM BUCKLES Sold in pairs.

- FS032-2** Width: 1" (25mm)
Ø Hole: 3/16" (5mm)
- FS033-2** Width: 1-1/2" (38mm)
Ø Hole: 1/4" (6mm)
- FS034-2** Width: 2" (50mm)
Ø Hole: 1/4" (6mm)



MOUNTING CLAMP

- HW404** Tube size: Ø3/4" (19mm) X Ø7/8" (22mm)
- HW404-2** Pair
- HW406** Tube size: Ø3/4" (19mm) X Ø1" (25mm)
- HW406-2** Pair

CLICKIT® AND FIXIT™ SEAT LATCH AND CATCH



CLICKIT® SEAT-LATCH

Spring-loaded latch release; secures solid seating drop hooks to wheelchair tubing.
Fits Tube Size: Ø1" (25mm) and Ø7/8" (22mm) with shim.

- HW600-1** Quantity: 1
- HW600-4** Quantity: 4
- HW600-10** Quantity: 10
- HW600-50** Quantity: 50
- HW600-100** Quantity: 100



FIXIT™ SEAT-CATCH

Secures solid seating drop hooks to wheelchair tubing.
Fits Tube Size: Ø1" (25mm) and Ø7/8" (22mm) with shim.

- HW601-1** Quantity: 1
- HW601-10** Quantity: 10
- HW601-50** Quantity: 50
- HW601-100** Quantity: 100



SEAT TUBE MOUNTING BRACKET KIT

These heat-treated, steel brackets allow for the quick installation and adjustment of hip belts directly to solid seats and rigid chair upholstery strips.

- HW312-KIT** Includes:
 2 Mounting Brackets
 4 10-32 x 3/4" Phillips Truss Screws

BAND CLAMP

This stainless steel, high-strength clamp fits into the tightest gaps, yet ensures a rock-solid hold of belts, harnesses and seating systems. Resistant to corrosion, it is ideal for use in showers and aquatherapy.

- HW320-22-2** Ø22mm (7/8")
- HW320-25-2** Ø25mm (1")
- HW320-27X33-2** Ø27mm X 33mm Oval (5/64" X 5/16")
- HW320-29-2** Ø29mm (1-5/32")
- HW320-30-2** Ø30mm (1-3/16")
- HW320-32-2** Ø32mm (1-1/4")

Above listed Band Clamp SKUs are also available in bulk packs of 10 (-10) and 50 (-50).




Band clamp has been tested to standards as a part of the complete postural support device mounting system.




BELT MOUNTING EXTENSION KIT

The Belt Mounting Extension Kit was created to give added length to the original Bodypoint® Hip Belts and the Evoflex® Pelvic Stabilizer while keeping all the same great features and fit.


- HKA12-1** Includes:
 4 Extension Straps
 2 M6 x 12 Button Head Socket Cap Screws
 2 M6 Nylock Nuts



Straps align to fit in tight spaces.



Straps spread apart to span obstructions.



Straps can be bent to attach directly to solid seat bottoms.

WARRANTY + TRADEMARK INFORMATION

Permobil Australia | Unit 3A/6 Boundary Rd | Northmead, NSW 2152 | www.permobil.com.au

LIMITED LIFETIME WARRANTY

Most products manufactured by Bodypoint carry a limited-lifetime warranty against defects in workmanship and materials arising under normal use by the original consumer. Product specifications are subject to change. For current specifications and warranty details, contact your local Independent Distributor or see our website at www.bodypoint.com.

LIMITATIONS ON WARRANTIES

The warranties described herein do not extend to products that have been subject to misuse, accident or tampering, nor to products that have been modified in any way other than that specifically described in the product instructions. Additional limitations may apply to distributed products not manufactured by Bodypoint – see specific manufacturer's warranty information.

Purchaser's sole remedy for defective products shall be repair, replacement or account adjustment at the discretion of Bodypoint. In all circumstances Bodypoint reserves the right to substitute an equivalent product or refund the purchase price of the original product.

The foregoing warranties are in lieu of all other warranties expressed or implied, including without limitation, implied warranties of merchantability and fitness for a particular purpose. In no event shall Bodypoint be liable for any special, indirect, incidental or consequential losses or damages resulting from the use of the product, including without limitation, loss of anticipated profits, even if Bodypoint shall have been advised of the possibility of such potential loss or damage. Bodypoint shall not be liable for any loss or damage caused by delay in furnishing product.

Bodypoint does not authorize any distributor, dealer, agent, sales person or any third party to change the terms or coverage of this warranty, and any such change shall be null and void. This warranty is non-transferable. This warranty will be governed in accordance with the laws of the State of Washington. The customer consents that jurisdiction under Venue for any action shall be in King County, Washington. This warranty information is reprinted here for information purposes only and is subject to change without notice.

TRADEMARK INFORMATION

Aeromesh is a registered trademark of Bodypoint, Inc.

Ankle Huggers is a registered trademark of Bodypoint, Inc.

Bodypoint is a registered trademark of Bodypoint, Inc.

® is a registered trademark of Bodypoint, Inc.

Clickit is a registered trademark of Bodypoint, Inc.

Cinch-Mount is a trademark of Bodypoint, Inc.

Evoflex is a registered trademark of Bodypoint, Inc.

Fixit is a trademark of Bodypoint, Inc.

Flat-Mount is a trademark of Bodypoint, Inc.

FrameSaver is a trademark of Bodypoint, Inc. |

Hip Grip is a registered trademark of Bodypoint, Inc.

Laminar Pad is a trademark of Bodypoint, Inc.

Monoflex is a trademark of Bodypoint, Inc..

PivotFit is a trademark of Bodypoint, Inc.

Rehab Latch is a trademark of Bodypoint, Inc.

Stayflex is a trademark of Bodypoint, Inc.

Taperlock is a trademark of Bodypoint, Inc.

MARK IV, V and VI are trademarks of Invacare Corporation.

Qtronix is a trademark of Sunrise Medical

Quickie is a registered trademark of Sunrise Medical

PGDT Pilot+ is a registered trademark of PG Drives Technology, Ltd. (formerly, Penny and Giles)

© Bodypoint® 2014