


DYNA-FORM™ Jupiter

High Risk Category



The Jupiter mattress can be used with patients considered to be at “High Risk” of pressure ulcer development.

The Jupiter design uses modern castellated foam on one side which provides excellent patient comfort and aids pressure ulcer management.

This mattress does not require turning which aids moving and handling procedures. It also works effectively with

any electric profiling bed frame.

The high quality CME foam and fabrics used ensure this modern designed mattress is durable and provides a solution with enhanced infection control.



Cross-section of castellations
(tear duct cut)

Castellated foam - tear duct cut

Smaller cells for heels - larger cells sacrally

Single side cut surface - works without turning programme

Multi-stretch, vapour permeable fabric - reduces shear and friction forces

Cover impregnated with anti-fungal agent

High frequency welded cover - prevents fluid ingress

Fully replaceable parts

	Width	Length	Height	Product Weight	Warranty	S.W.L.
Single	85.5cm	198cm	15cm	10kg	Cover - 3 years Foam - 3 years	165kg
King Single	107cm	198cm	15cm	12kg	Cover - 3 years Foam - 3 years	165kg