

DYNA-FORM™ Mercury Advance Bari

Very High Risk Category

The Mercury Advance Bari is specifically designed for patients at "Very High Risk" of pressure ulcer development. It combines the benefits of modern foam technology with the facility to "step up" to a dynamic mattress when clinically required. This makes the Mercurv Advance the ideal system to use Foam within in a patient's home or palliative care individual air cell. Cells can be replaced environment as it helps reduce logistic and separately. decontamination costs.

It has a higher maximum weight capacity of 380kg to meet the needs of heavier clients. All component parts are interchangeable and replaceable, maximising product life and reducing environmental impact.

The outer cover comprises a high-frequency welded, multi stretch and vapour permeable fabric which satisfies the strictest infection control policies.



CME foam castellations (Cross-Section)

Fixed head section - head | shoulder stability Static mode - rapid transfer or cleaning Pressure selection mode - patient comfort 2-in-1 mattress - minimises logistics and decontamination Easy use power unit - automatic pressure adjustment

Foam and dynamic mattress combination

LED Alert - flashing indicators for low cell pressures or failure

CPR connector valve - rapid deflate within 10 seconds

	Width	Length	Height	Product Weight	Power Unit Weight	Warranty	S.W.L.
Bari	120cm	198cm	15cm	25kg	1.7kg	Cover - 4 years Foam - 8 years Pump Cells - 2 years (excludes pump bellows - 12 months)	380kg