

PRODUCT DATASHEET

Disposable Jug

Brand	: CURAS [®]
Category	: Nursing
Range	: Hygiene Pulp [™]
Product Group	: Human Waste Collection
Product REF	: 08034
Product Description	: Disposable Jug 1L, WardPak [™]



PRODUCT DESCRIPTION

Disposable jug made of pulp material with water resistance properties. 1L capacity. Suitable for disposal in a macerator. Packed in WardPak[™] for easy dispensing. Non-sterile.

INTENDED USE

Disposable jug for draining urine from urine drainage bag/system. It is also suitable for holding other liquid contents. Single use.

FEATURES & BENEFITS

- An environmental friendly product that is made of 100% recycled newsprint paper.
- Disposable for reduced risk of cross-contamination and Hospital-Acquired Infections.
- Water resistant, for a minimum period of 4 hours (applicable to jug only).
- Integrated graduation marking for indication of approximate volume.
- Reinforced rim provides extra strength and supports weight of contents when full.
- Ergonomic design allows easy holding and handling.
- Rounded edges on the base prevent cracking and increase wet strength.
- Reinforced base provides extra stability.
- Double-press process gives a smooth surface and better finishing on product.
- Stackable to reduce storage space.
- Inner bag, WardPak[™] allows easy dispensing at ward level and protects product from dust and contamination.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	Packing
08034	1L	25.8 x 18.8 x 9.98 cm	40g ± 3 g	WardPak™

Composition

- Jug is made of 100% recycled newsprint paper with water-based fluid resistant compound.

Country of Origin

- China

Sterilization

- Non-sterile.
- Can be sterilized by EtO and Gamma radiation.

Shelf Life

- 8 years

Selected Standards & Guidelines

- ISO 9001:2015
- MDR 2017/745, Class I
- ISO 14971:2012
- ISO 20417:2021
- ISO 15223-1:2012
- MEDDEV 2.12-1 rev 6 (2009)
- EN 15223-1:2021
- PAS 29:1999

Product Symbols



Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose of used product into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

PACKAGING INFO

REF	GTIN-14	WardPak™	Carton	Gross Weight
08034	05051400080340	100 pcs / bag	Jug: 200 pcs / carton	~9 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	PE	--	✓	(✓)		✓	
Outer Carton	Corrugated Carton	--	(✓)	(✓)		✓	
Transport Pallet	Wood	--				(✓)	✓

✓ = Highly recommended

(✓) = Recommended

ADDITIONAL INFO

- Curas[®] Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com