

PRODUCT DATASHEET

Disposable Ear Irrigation Cup

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product Group	: Human Waste Collection
Product REF	: 08036
Product Description	: Disposable Ear Irrigation Cup (Noots Tank), 300ml, WardPak™



PRODUCT DESCRIPTION

Disposable ear irrigation cup (noots tank), also known as ear flush cup, is made of pulp material with water resistance properties. 300ml capacity. Suitable for disposal in a macerator. Packed in WardPak™ for easy dispensing. Non-sterile.

INTENDED USE

Disposable cup for collection of fluids and waste during the ear cleansing procedure. Single-use only.

FEATURES & BENEFITS

- An environmental friendly product that is made of 100% recycled newsprint paper.
- Safe and convenient option to collect ear flushing fluids.
- Water resistant, for a minimum period of 4 hours
- Depth of more than 14 cm for high fluid capacity minimizes spills and splashes.
- Shaped to fit under the earlobe for an easy access flushing procedure.
- Cup shaped for easier handling.
- Hygienic due to one patient usage, ensuring a new cup for each user.
- Does not require washing and disinfecting, which is more economical than non-disposable.
- Can be conveniently disposed with waste fluid into the macerator.
- Stackable for reduced storage space.
- Inner bag, WardPak™ protects product from dust and contamination.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (Ø x H)	Average Weight/Pc	Packing
08036	300ml	Ø 8.5 x 14.5 cm	20 g	WardPak™

Composition

- Made of 100% recycled newsprint paper.
- Water-based fluid resistant compound.

Country of Origin

- China

Sterilization

- Non-sterile.

Shelf Life

- 8 years

Selected Standards & Guidelines

- ISO 9001:2015
- MDR 2017/745, Class 1
- ISO 13485:2016
- ISO 14971:2019
- ISO 20417:2021
- ISO 15223-1:2021
- PAS 29:1999

Recommended Storage Condition

- Store at dry, clean, and well-ventilated area at room temperatures between 10° to 30°C (50°-86°F).
- Keep away from moist and store at dry place with humidity level between 40-60%RH.
- Avoid exposure to extreme temperatures. Keep away from moisture, heat and sources of ignition.
- Do not expose the product to radiation, sunlight and UV light, as it will change the properties of the product.

Proposed Disposal Method

- If regular waste, dispose of used product into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

Product Symbols



PACKAGING INFO

REF	EAN-13 Barcode	WardPak™	Carton	Gross Weight
08036	5051400080364	40 pcs / bag	160 pcs / carton	~3.09 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√= Highly recommended (√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com