

## PRODUCT DATASHEET

# Disposable Urine Bottle, Female

800ml

Brand : CURAS® Category : Nursing

Range : Hygiene Pulp™

Product Group : Disposable Urine Bottles

REF : 08052



#### PRODUCT DESCRIPTION

Disposable female urine bottle made of hygienic pulp material with water resistance properties. Suitable for disposal in a macerator. 800ml capacity.

#### **INTENDED USE**

Disposable urine bottle is used by immobile female patients for collection of urine and draining of urine bags. Single use.

#### **FEATURES & BENEFITS**

- Disposable for reduced risk of cross-contamination and Hospital-Acquired Infections (HAIs).
- Water resistant, for a minimum period of 4 hours.
- Ergonomic design for easy handling and convenient to use.
- Product is steady with flat bottom.
- Fluid Absorber, SAP (REF 08100) can be inserted prior to use for improved waste handling.

#### **PRODUCT SPECIFICATIONS**

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	Packing	
08052	800ml	25 x 11 x 10.8 cm	60 g ± 2 g	Carton	

#### Composition

- Made of 100% recycled newsprint paper.
- Water-based fluid resistant compound.

#### **Country of Origin**

China

## Sterilization

- Non-Sterile.
- Can be sterilized by EtO and Gamma radiation.

## **Shelf Life**

8 years



#### **Selected Standards & Guidelines**

- MDD 93/42/EEC (CE Marking, Class I), amended by Directive 2007/47/EC.
- ISO 9001:2008
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at 35±3°C).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

## **Recommended Storage Condition**

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

## **Proposed Disposal Method**

- If regular waste, dispose of used product into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

## **Product Symbols**













curas.com/Symbols

#### **PACKAGING INFO**

REF	EAN-13 Barcode	Inner Bag	Carton	Gross Weight	
08052	5051400080524	-	21 pcs / carton	~1.8 kg	

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
Outer Carton	Corrugated Carton		(√)	(√)		√	
Transport Pallet	Wood					(√)	√

 $<sup>\</sup>sqrt{\ }$  = Highly recommended

## **ADDITIONAL INFO**

- Curas® Hygiene Pulp™ products have been tested in several Haigh macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com

 $<sup>(\</sup>sqrt{\ })$  = Recommended