

PRODUCT DATASHEET

Disposable Female Urine Bottle

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product REF	: 08052
Product Description	: Disposable Female Urine Bottle, 500ml, WardPak™



PRODUCT DESCRIPTION

The Disposable Female Urine Bottle is made of a pulp material that has water resistant properties. Its volume capacity is 500ml. It is suitable for disposing into a macerator.

INTENDED USE

The Disposable Female Urine Bottle is used by immobile female patients for the collection of urine. It is for single-use only.

FEATURES & BENEFITS

- Reduces the risk of cross-contamination and Hospital-Acquired Infections (HAI) as this is a single-use disposable product.
- Water resistance ensures safe use with liquid for at least 4 hours.
- Ergonomic design which gives a good fit to a woman's external anatomy.
- Comes with a contoured handle for ease of gripping the bottle and handling.
- Convenient to use and provides privacy to female patients. Patients are able to have toileting independence without the help of a nurse.
- Can be used with a SAP Pad by Curas™ (REF 09602) to prevent accidental spills, control odour and improve waste handling.
- This product has both smooth surface and finishing that is clean to use.
- This product is steady due to its flat bottom.
- Able to fit into common bedside holders.
- WardPak™ protects product from dust and contamination. The perforated opening on the dispensing carton also offers an easy dispensing option.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H) cm	Average Weight/Pc	Packing
08052	500ml	22 x 12 x 9	50 g ± 3 g	Carton

Composition

- Made of recycled newsprint paper with water-based fluid-resistant compound.

Country of Origin

- China

Sterilization

- Non-sterile.
- Can be sterilized by EtO (gas) and Gamma (radiation)

Shelf Life

- 8 years

Selected Standards & Guidelines

- ISO 9001:2015
- MDR 2017/745, Class I
- ISO 14971:2012
- BS ISO 20417:2021
- ISO 15223-1:2016
- PAS 29:1999

Recommended Storage Condition

- Store at room temperature (15°C-25°C / 59°F-77°F) and avoid exposure to elevated temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method



- If regular waste, dispose used product into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or dispose in a landfill.
- If content is contaminated or contagious, treat as clinical waste and dispose in accordance with the clinical waste protocol.

Product Symbols



PACKAGING INFO

REF	GTIN	Inner bag	Carton	Gross Weight
08052	05051400080524	-	24 pcs / carton	~ 2 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Symbol	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	White Corrugated Carton		(√)	(√)		√	
Transport Pallet	Wood					(√)	√

√ = Highly recommended (√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com