

Rоно[®] Quadtro Select[®] Series

Very High Risk Category



Roho Dry Floatation cushions constantly adjust and conform with an individual's body and support the shifts and movements made throughout the day. The revolutionary ISOFLO[®] Memory Control offers customized shape-fitting capabilities while you're seated with just the push of a button. It allows for maximum skin protection without compromising on stability, positioning or convenience.

The individual's posture can be adjusted to the desired sitting position by regulating air in four quadrants and then the positioning settings saved by pushing the knob to the closed position.

The Quadtro series cushions are available in different profiles to accommodate varying degrees of pelvic asymmetry and pressure care requirements.

Made of flame resistant neoprene rubber Includes a two-way stretch cover, hand inflation pump and repair kit

Neoprene rubber cushion can be patched or repaired to extend use

Shapes to the individual, addressing their specific seating needs

Adds little weight to the full mobility system

Cushion is easily cleaned with soap and water.

Cover is machine washable.



ISOFLO Control **Open** - Air can flow freely throughout entire cushion ISOFLO Control **Closed** - air becomes locked into four separate quadrants

	Width	Depth	Height	Product Weight	Risk Category	Warranty	S.W.L.
Low Profile	28 - 62.5cm	30 - 56cm	5.5cm	1.1kg	Very High	2 years	No limit but cushion must be properly sized to individual
High Profile	32.5 - 67	30 - 51.5cm	10.5cm	1.5kg	Very High	2 years	

WA (08) 9456 1661

VIC (03) 9580 1055

SA (08) 7225 8720 www.hospequip.com.au

QLD (07) 3448 9243

NSW 1300 467 737