
VICAIR® ACTIVE 02

FOR THE ACTIVE WHEELCHAIR USER









PRESSURE ULCER TREATMENT

Excellent lateral and frontal stability

Perfect fit for the user who transfers frequently

Superior ischial offloading and protection

The Active O2 offers the unique combination of ergonomically shaped foam and Vicair SmartCells. The four rear SmartCell compartments provide a high level of comfort and skin protection, and offload the sitting bones and tailbone area. The foam front makes transfers easy and feels steady and improves wheelchair handling. The foam inlay is protected by an easy-to-clean fluidresistant liner. The rest of the cushion is 100% machine washable and breathable for optimal hygiene and perfect temperature & moist regulation. The Active O2 has an easy accessible storage pouch for precious belongings (keys, phone, wallet). In short: everything you need to maintain your active lifestyle!

PRODUCT PERFORMANCE		
Skin protection / Anti decubitus	***	
Pressure redistribution	***	
Shear reduction	***	
Micro climate control		
Impact & vibration damping	***	

Positioning	***
Lateral stability	AAAA
Frontal stability	AAAA
Sliding prevention	AAAA
Allows manual adjustments	**
Allows easy transfers	***

Maintenance & Safety	***
Durability	AAAA
Reliability	AAAA
Functionality when misused	
Low maintenance level	
Level of hygiene	***

Functionality & Application	***
Comfort	AAAA
Low cushion weight	
Available size range	**
Maximum user weight	**
Low seat to floor height	AAAA
Auto adjustments	

USED WITH

Spinal cord injuries, active users with a high transfer frequency

TECHNICAL SPECIFICATIONS	
Height	6 and 9 cm
Average weight	700 grams
Max. user weight	200 kg
Machine washable	60°C (without foam)

CUSHION SIZES	
Standard	(Sizes >53 cm height 9 only)
On request	

