VICAIR® MULTIFUNCTIONAL 02

COMFORT AND SKIN PROTECTION IN A TILT-IN-SPACE WHEELCHAIR

Longer compartments to fit into tilt-in-space chairs | follow the body's contours for equal pressure distribution | support IT-coccyx area | managing pelvic asymmetry









PRESSURE ULCER TREATMENT

Ideal for use in tilt-in-space chair with adjustable seat depth

Vicair Adjuster O2 design with deeper rear end compartments

Optimal adjustment to the body's contours for equal pressure distribution

The Vicair Multifunctional O2 seat cushion is developed to use in a tilt-in-space wheelchair with adjustable seat depth. It provides a high level of comfort, optimal skin protection and stability. The design allows total immersion to redistribute pressure equally, and it keeps you comfortably positioned and supported. 100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation. A safe and reliable solution, especially when used together with the Vicair

PRODUCT PERFORMANCE	
Skin protection / Anti decubitu	S AAAA
Pressure redistribution	***
Shear force reduction	AAAA
Microclimate control	AAAA
Impact & vibration damping	AAAA
Positioning	***
Lateral stability	
Frontal stability	AAAA
Sliding reduction	***
Allows manual adjustments	***
Allows easy transfers	
Maintenance & Safety	***
Durability	***
Reliability	***
Functionality when misused	
Low maintenance level	***
Level of hygiene	***
Functionality & Application	***

Comfort

Low cushion weight

Available size range

Auto adjustments

Maximum user weight

Low seat to floor height

USED WITH

Neuromuscular diseases, Rheumatism, Quadriplegia, pain relief, in a geriatric environment

TECHNICAL SPECIFICATIONS	
Height	8 cm
Average weight	850 grams
Max. user weight	250 kg
Machine washable	60°C max.

CUSHION SIZES	
cm	inches
37x52	14¾x20¾
40x52	16x20¾
45x52	18x20¾
50x52	20x20¾

Sizes in width x depth