


MATRIX® Absolute Cushion



Simple, highly functional solution for individuals with mild to moderate seating needs

Highly resilient molded polyurethane foam

Gentle contour for comfort and centring

Water repellent outer cover protects foam

For individuals still able to do own weight shift to provide pressure relief

Shaped cushion bottom facilitates use on taut sling upholstery

Hook and loop fasteners keep cushion securely in place

Seat Rigidizer

Lightweight plastic seat insert

For use on sling seat wheelchairs

Eliminates hammock effect of sling seat upholstery

Fits between seat rails for level support


Does not change weight limitation of wheelchair

Securely fitted to underside of cushion

foam using double sided tape

Fits inside cushion cover

Available in various widths and depths to accommodate most cushion sizes



	Width	Depth	Height	Product Weight	S.W.L.
Absolute	40.5 45.5 51cm	40.5 45.5cm	Adductor - 8cm Abductor - 5.5cm	0.9kg	113kg
Rigidizer	40 45 50	40 45cm	-	-	-