

Freedom Turning Bed



The Freedom Turning Bed is unique in its design, with features to accommodate the needs of immobile individuals who are at risk of developing pressure injuries or who require a solution to treat existing pressure ulcers or complications caused by chronic immobility.

It is automated, programmable and has the option to be fully voice controlled.

The need for manual turning is eliminated allowing an increase in home-based care where previously institutionalization would have been necessary. In addition, both the patient and the carer experience improved outcomes from subsequent undisturbed sleep.

Mimics natural nocturnal motion by repositioning to re-establish blood flow to compressed tissue

Lower side of mattress turns up to form a “wing” which supports along entire body length from ankle to shoulder

Turning occurs along spinal axis
- minimal feeling of rotation

Air powered head and leg elevations

Torso position can be raised to 30° - rotation still possible with raised torso


Automatic or manual lateral rotation



30° both ways



3-sectioned wound care mattress



Trendelenburg and Reverse Trendelenburg

	Overall Width Mattress Platform	Overall Length Mattress Platform	Height range	Rotation angle	Torso raise angle	Trendelenburg tilt	S.W.L.
Single	110cm 97cm	262cm 218cm	56 - 86cm	30° - 0° - 30°	0° - 30°	0° - 30°	159kg (272kg - Bari)