

DYNA-FORM SMARTRESPONSE MATTRESS

Dyna-Form® SMARTresponse® is a next generation hybrid support surface that is uniquely able to support clinical choice and timely intervention on an individualised basis by automatically or remotely stepping-up patient care when a period of patient non-movement has been detected. With the focus on clinical choice, SMARTresponse® empowers healthcare providers to achieve earlier intervention in pressure ulcer prevention.

Hybrid support surfaces have been clinically proven to support the prevention and treatment of pressure ulcers, representing a step-change in care delivery that provides a more cost-effective solution and promotes best use of resources. With the addition of the SMARTresponse® control unit, DHG have transformed the hybrid support surface, providing a unique smart-monitoring mattress that can also provide intervention at the point of increased risk.

SMARTresponse® offers the following unique smart-technology features:

- **Auto Start** – Selecting Auto Start allows the system to automatically ‘step up’ care when no patient movement has been detected for a period of time.
- **Alert Me** – Alert Me sends a message directly to the care provider when no patient movement has been detected for a period of time. This provides the carer with the ability to remotely ‘step up’ care from anywhere.
- **Bed Exit** – Bed Exit can be activated at any point to detect bed exit events. In the event that a bed exit is detected, an alert is sent directly to the care provider.



Auto Start



Alert Me



Bed Exit



254KG SWL



VERY
HIGH RISK
CATEGORY



SUITABLE FOR
STAGE 4
MANAGEMENT



2 YR FOAM
4 YR COVER
WARRANTY



CONTACT US



SMARTresponse transforms the pressure ulcer prevention care delivery model, providing the next-generation hybrid mattress.



Hybrid support surfaces are clinically proven to support the prevention and treatment of pressure ulcers. Representing a step change in care delivery, hybrid solutions provide a more cost-effective solution, promoting best use of resources when compared to the traditional use of standard high specification foam and alternating pressure air mattresses.

With the addition of the **SMART**response power unit, DHG have transformed the hybrid support surface, enabling the immediate 'step up' of patient care automatically or remotely when a period of no patient movement has been detected:



CONTINUALLY OPTIMISED
PRESSURE REDISTRIBUTION

NON-MOVEMENT MONITORING
FOR EARLIEST INTERVENTION

EFFECTIVE PRESSURE RELIEF

SPECIFICATIONS

SINGLE SKU: 6DHG-MAS321002096 | KING SINGLE SKU: 6DHG-MAS3210001

Model	Width	Length	Thickness	Product Weight	Warranty	SWL
Single	88cm	198cm	15cm	13.4kg Mattress (1.7kg control unit)	2 Years	254kg
King Single	107cm	198cm	15cm			

CONTACT US