



EQUAGEL™ General and Protector Cushion



EquaGel cushions are designed with two layers of rubbery dry-polymer gel. It is not a liquid. It works on a principle referred to as column buckling. Each gel wall can only support a marginal load. When that threshold is exceeded, it buckles under and passes the weight to the surrounding walls, thereby effectively distributing the pressure by maximising the surface area.

Available in a range of sizes, in either a 5cm General or 6.3cm Protector option. The top comfort layer comprises thin wall gel in a tight grid and the bottom layer is thicker wall gel with a wider grid for deep posture support. The contoured pelvic area provides additional pressure distribution while maintaining gel below the pressure points.

- Made of 100% EquaGel
- Equagel is a dry-polymer gel
- 5cm and 6.3cm height - in two stages
- Contoured pelvic area distributes weight evenly
- Retains buckling gel under the pressure points
- High level of heat and vapour permeability
- Non-allergenic - waterproof, easily washable
- Fluid resistant, vapour permeable cover
- Machine washable cover



Tight grid top comfort layer and wide grid deep posture support layer

| | Width | Depth | Height | Risk | S.W.L. |
|-----------|-----------------------|---------------|--------|------------|--------|
| General | 40.5 45.5 51 56 | 40.5 45.5cm | 5cm | Low | 125kg |
| Protector | 61 66 71 76 cm | | 6.3cm | Med High | 147kg |