


### DYNA-FORM MERCURY ADVANCE MATTRESS



254KG SWL



VERY HIGH RISK



SUITABLE FOR UP TO STAGE 4 MANAGEMENT



8 YEAR FOAM WARRANTY



**CONTACT US**

#### FEATURES

#### BENEFITS

Combined foam and dynamic mattress technology

Provides the ability to 'step up' and 'step down' care to that of the users clinical requirements

Fixed head section

Maintaining head and shoulder stability

Static & Pressure Selection mode

Static for rapid transfer or cleaning procedures, Pressure for providing ultimate patient comfort, without affecting care procedures

2-in-1 mattress capability

Minimises demands on logistics, decontamination and repair costs

Clinically proven foam mattress

Over 25 clinical studies on the Mercury Advance

Audible-visual warning system

Flashing LED and audible warning system alerts the carer/individual to low cell pressure or system failure

Easy deflate CPR connector valve

Allows rapid deflation within 10 seconds

The Mercury Advance is specifically designed for patients at “Very High Risk” of pressure ulcer development. It combines the benefits of modern foam technology with the facility to “step up” to a dynamic mattress when clinically required. This makes the Mercury Advance the ideal system to use in a patient’s home or palliative care environment as it helps reduce logistic and decontamination costs. It has a maximum weight capacity of 254kg to accommodate heavier patients.

All parts are interchangeable and replaceable, maximising product life and reducing environmental impact. The outer cover comprises a high-frequency welded, multi-stretch and vapour-permeable fabric which satisfies the strictest infection control policies. The pump has flashing LED alert indicators to alert carers to low cell pressures or failures. The CPR connector valve allows for rapid deflation (within 10 seconds) of the air cells in the event CPR is necessary. Pump sold separately.

### KEY FEATURES

- Foam and dynamic mattress combination
- Fixed head section - head | shoulder stability
- Static mode - rapid transfer or cleaning
- Pressure selection mode - patient comfort
- 2-in-1 mattress - minimises logistics and decontamination
- Easy use power unit - automatic pressure adjustment
- LED Alert - flashing indicators for low cell pressures or failure
- CPR connector valve - rapid deflate within 10 seconds



### SPECIFICATIONS

Model	Width	Length	Thickness	Product Weight	Foam Density	SWL	Warranty
Single	88cm	198cm	15cm	13.3kg	Foam Insert - 38-40kg/m3 Nominal Hardness - 120N	254kg	8 Year Foam 4 Year Cover Pump & Cells 2 years
King Single	107cm	198cm	15cm	15kg	Blue CME Foam - 38-40kg/m3 Nominal Hardness - 200N	254kg	Pump bellows 12 months

**CONTACT US**