

JAY[®] Ion Cushion



The Ion cushion features a dual-layered foam base. The top layer consists of Visco memory foam which conforms to the individual and distributes their weight over a larger area to reduce pressure and maximise comfort. To optimise immersion, the foam's firmness is determined by the cushion width, which directly corresponds to the client's weight.

The foam base is contoured to help redistribute the load and also contribute to correct positioning.

The dual cover system increases the longevity of the cushion by protecting it from fluids through the water resistant inner cover, which is easy to wipe off and clean. The stretch outer cover is machine-washable and utilizes silver-impregnated X-static[®] fibre which is naturally anti-microbial.

- Dual layer foam base - varying firmness
- Visco memory foam top layer
- Pre-contoured base helps redistribute load
- Dual cover system
- Reticulated foam in outer cover promotes air movement
- Inner cover is water resistant and easy to wipe off
- Aquaguard[™] zipper resists moisture
- Stretchable outer cover
- Additional outer covers available separately



Dual cover system helps protect the cushion.

The quick dry, machine washable outer cover uses X-static[®] stretch fabric which is impregnated with silver fibres which naturally inhibit bacteria growth

Must be sized correctly to individual	Width	Depth	Height	Product Weight	S.W.L.
	35.5 - 61cm	35.5 - 51cm	10cm	1.5kg (varies by size)	135kg (W 35.5 - 53cm) 225kg (W 56 - 61cm)